

USER MANUAL

InviMag[®] Plant DNA Mini Kit/ KF96

for semi-automated purification of genomic DNA from up to 100 mg plant material and food of vegetable origin with magnetic beads on the on KingFisher™ 96 and KingFisher™ Flex, Thermo Fisher Scientific

Instruction for InviMag® Plant DNA Mini Kit/ KF96

The **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** is the ideal tool for isolation and purification of DNA from up to 100 mg plant material or food from plant origin with the patented InviMag® technology using a KF96 / KFflex96 instrument.

The kit is neither validated for the isolation of genomic DNA from cultured or isolated cells, from blood, stool samples, swabs, dried blood stains or cell free body fluids, like cerebrospinal fluid, synovial fluid and urine nor from bacteria, fungi, parasites or purification of RNA.

Trademarks: InviMag®, Invisorb®. Registered marks, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

The Invisorb® technology is covered by patents and patent applications: US 6,110,363, US 6,043,354, US 6,037,465, EP 0880535, WO 9728171, WO 9534569, EP 0765335, DE 19506887, DE 10041825.2, WO 0034463.

InviMag® and Invisorb® are registered trademarks of Invitex Molecular GmbH.

The PCR process is covered by US Patents 4,683,195, and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG.

© 2020 Invitex Molecular, all rights reserved.

Table of content

Kit contents of InviMag® Plant DNA Mini Kit/ KF96 w/o plastic	3
Symbols	4
Storage	4
Intended use	5
Product use limitation	5
Safety information	6
Product characteristic of the InviMag® Plant DNA Mini Kit/ KF96 w/o plastic	7
Sampling and storage of starting material	8
Principle and Procedure	8
Yield and quality of genomic DNA	8
Important points before starting a protocol	9
Preparing reagents and buffers	9
Reagents and equipment to be supplied by user	9
Scheme of the InviMag® Plant DNA Mini Kit/ KF96 w/o plastic	10
Lysis Procedures	11
<i>Protocol: Isolation of genomic DNA from up to 100 mg of plant material</i>	11
Starting a Run on a KF96 / KFflex96 instrument	11
For self-programming of the KF96 / KFflex96 instrument	13
Troubleshooting	15
Appendix	16
General notes on handling DNA	17
Ordering information	18
Ordering information (KingFisher™ 96 and consumables)	18

Kit contents of InviMag® Plant DNA Mini Kit/ KF96 w/o plastic

	5 x 96 extractions
Catalogue No.	7437300250
Lysis Buffer P	210 ml
Proteinase K working solution	10 x 1.1 ml working solution
Binding Buffer A	36 ml (final volume 120 ml)
SNAP Solution	10.5 ml
Wash Buffer I	3 x 80 ml (final volume 3 x 160 ml)
Wash Buffer II	4 x 60 ml (final volume 4 x 200 ml)
Elution Buffer	60 ml
Sealing Foils	10
Manual	1
Prefilter Plate I	5
Initial steps	<p>Resuspend each tube Proteinase K in 1.1 ml RNase free water, mix thoroughly until completely</p> <p>Add 80 ml of 96-100% ethanol to each bottle Wash Buffer I</p> <p>Add 140 ml of 96-100% ethanol to each bottle Wash Buffer II, mix thoroughly and always keep the bottles firmly closed</p> <p>Add 84 ml 99.7% Isopropanol to the Binding Buffer A. Mix by intensive shaking by inverting for 1 min. Shortly before use mix by inverting several times.</p>
Plastic to be supplied by user (see order information)	
2.0 ml Deep Well Plate	20
KF 96 Tip Comb for DW magnets	5
200 µl Elution Plate*	10

*) Elution and Tip Plate are identically. Use one provided Elution Plate as a Tip Plate.

Symbols



Manufacturer



Lot number

Attention: Do not combine components of different kits, unless the lot numbers are identical!



Catalogue number



Expiry date



Consult operating instructions



Temperature limitation



Do not reuse



Humidity limitation

Storage

All buffers and kit contents of the **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic**, except **dissolved Proteinase K** should be stored at room temperature and are stable for at least 12 months.

Room temperature (RT) is defined as range from 15-30°C.

Before every use, make sure that all components have room temperature. If there are any precipitates within the provided solutions solve these precipitates by warming carefully (up to 30°C).

Proteinase K: Dissolved Proteinase K must be stored at 2 - 8 °C for up to two months. For longer storage -20 °C is recommended, freeze-thaw once only. Therefore, the dissolved Proteinase K is stable as indicated on the kit package.

Wash Buffers charged with ethanol should be appropriately sealed and stored at room temperature.

Binding Buffer charged with isopropanol should be appropriately sealed and stored at room temperature.

Quality control and product warranty

Invitek Molecular warrants the correct function of the **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** for applications as described in this manual. Purchaser must determine the suitability of the product for its particular use. Should any product fail to perform the applications as described in the manual, Invitek Molecular will check the lot and if Invitek Molecular investigates a problem in the lot, Invitek Molecular will replace the product free of charge.

Invitek Molecular reserves the right to change, alter, or modify any product to enhance its performance and design at any time. In accordance with Invitek Molecular's EN ISO 13485 certified Quality Management System the performance of all components of the **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** have been tested separately against predetermined specifications routinely on lot-to-lot to ensure consistent product quality.

If you have any questions or problems regarding any aspects of **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** or other Invitek Molecular products, please do not hesitate to contact us. A copy of Invitek Molecular's terms and conditions can be obtained upon request or are presented at the Invitek Molecular webpage www.invitek-molecular.com.

For technical support or further information please contact:

from Germany +49-(0)30-9489-2901/ 2910

from abroad +49-(0)30-9489-2907

or contact your local distributor.

Intended use

The **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** is designed for semi-automated preparation of genomic DNA from up to 100 mg plant or food (plant origin) material using the patented **InviMag® technology** in combination with a KF96/ KFflex96 instrument.

The whole process is based on a patented bead technology, used for isolation of genomic DNA by binding the nucleic acid onto magnetic particles in absence of chaotropic buffer components.

For reproducible and high yields, appropriate sample storage is essential.

THE PRODUCT IS INTENDED FOR USE BY PROFESSIONALS ONLY, SUCH AS TECHNICIANS, PHYSICIANS AND BIOLOGISTS TRAINED IN MOLECULAR BIOLOGICAL TECHNIQUES. It is designed to be used with any downstream application employing enzymatic amplification or other enzymatic modifications followed by signal detection or amplification. Any diagnostic results generated by using the sample preparation procedure in conjunction with any downstream diagnostic assay should be interpreted with regard to other clinical or laboratory findings.

To minimize irregularities in your results, adequate controls for downstream applications should be used.

Product use limitation

The kit is neither suitable for the isolation of DNA from blood, serum or plasma, bacteria, fungi or viruses, nor for isolation and purification of RNA.

The included chemicals are only useable once.

Differing of starting material or flow trace may lead to inoperability; therefore, neither a warranty nor guarantee in this case will be given, neither implied nor express.

The user is responsible to validate the performance of the Invitek Molecular product for any particular use. Invitek Molecular does not provide for validation of performance characteristics of the product with respect to specific applications.

Invitek Molecular products may be used e.g. in clinical diagnostic laboratory systems under following conditions:

- If used in the US, based on the condition that the complete diagnostic system of the laboratory has been validated pursuant to CLIA' 88 regulations.
- For other countries based on the condition that the laboratory has been validated pursuant to equivalentents according to the respective legal basis.

All products sold by Invitek Molecular are subject to extensive quality control procedures (according to EN ISO 13485) and are warranted to perform as described herein. Any problems, incidents or defects shall be reported to Invitek Molecular immediately upon detection thereof.

The chemicals and the plastic parts are for laboratory use only; they must be stored in the laboratory and must not be used for purposes other than intended.

The product with its contents is unfit for consumption.

Safety information

When and while working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles!

Avoid skin contact! Adhere to the legal requirements for working with biological material!

For more information, please consult the appropriate material safety data sheets (MSDS). These are available online in convenient and compact PDF format at www.invitek-molecular.com for each Invitek Molecular product and its components. If buffer bottles are damaged or leaking, **WEAR GLOVES, AND PROTECTIVE GOGGLES** when discarding the bottles in order to avoid any injuries.

Invitek Molecular has not tested the liquid waste generated by the **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** procedures for residual infectious materials. Contamination of the liquid waste with residual infectious materials is highly unlikely, but cannot be excluded completely. Therefore, the liquid waste has to be considered infectious and be handled and discarded accordingly to local safety regulations.

European Community risk and safety phrases for the components of the **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** to which they apply, are listed below as follows:

Lysis Buffer P



Warning

H319-H412-P280-P305+P351+P338-P273

Proteinase K:



Danger

H315-319-334-335 P280-P305-P351-P338

Wash Buffer I



Warning

H302-412-P280-P305-P351-P338-P273-EUH032

H302: Harmful if swallowed.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H335: May cause respiratory irritation.

P273: Avoid release to the environment.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

EUH032: Contact with acids liberates very toxic gas.

Emergency medical information can be obtained 24 hours a day from infotrac:

outside of USA: 1 – 352 – 323 – 3500

inside of USA : 1 – 800 – 535 – 5053

Product characteristic of the InviMag® Plant DNA Mini Kit/ KF96 w/o plastic

Starting material	Yield	Time	Ratio
up to 100 mg plant material	5 - 25 µg; depends on the kind and amount of starting material	about 30 min (without lysis)	A ₂₆₀ :A ₂₈₀ 1.8-2.2

The **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** is designed for a semi-automated preparation of genomic DNA from up to 100 mg of plant material using magnetic beads and the KF96 or KFflex96 workstation. The isolation process is based on the patented **InviMag®** technology.

The DNA isolation process relies on the interaction of nucleic acids with coated magnetic particles at adapted buffer conditions. The KF96/KFflex96 instrument performs all purification steps of the DNA purification procedure automatically, except the plant sample preparation, lysis and initial loading of the system. Sample cross-contamination and reagent cross-over is effectively eliminated by the provided assay file.

The KingFisher® instrument uses magnetic rods to transport the DNA bound to magnetic particles through the various assay phases like binding, washing, drying and elution. The volume of buffers and other liquids required for DNA isolation is reduced to a minimum.

To achieve optimal lysis conditions and high yields, the plant samples are first mechanically disrupted followed by a lysis step in an optimized buffer system at elevated temperature. After lysis, a binding step is performed in which the DNA is bound to the magnetic particles followed by several washing steps before the pure DNA is finally eluted.

The purified and high quality DNA can be stored at -20°C for subsequent use and is ready-to-use for subsequent downstream applications like:

- PCR*
- Genotyping
- Restriction digestion

For further information please contact +49 (0) 30 9489 2901 or 2910 in Germany and +49 (0) 30 9489 2907 from foreign countries or ask your local distributor.

Sampling and storage of starting material

Harvested plant samples can be stored at room temperature for up to 2–3 hours. For short-term storage (up to one week), samples may be stored at 2-8°C. For long-term storage, we recommend freezing samples at -20°C or -80°C. Multiple thawing and freezing cycles before isolating the DNA should be avoided because this can lead to degraded DNA and reduced yields.

Invitex Molecular will be released of its responsibilities if other sample materials than described in the Intended Use are processed or if the sample preparation protocols are changed or modified.

Principle and Procedure

Lysis

Samples are disrupted by using a mixer mill, bead mill or by mortar and pestle in combination with liquid nitrogen. Afterwards, the disrupted material is lysed at denaturing non-chaotropic conditions at elevated temperatures in presence of **Lysis Buffer P** and **Proteinase K**. Because unlysed material has to be removed after lysis by filtration, the lysis procedure has to be performed externally. The filtration step is required, because in most cases, a centrifugation step will not be able to remove very small debris, especially if a mortar and pestle was used for disrupting.

Binding of the DNA

After addition of **Binding Buffer A** and **SNAP Solution** to the lysate, optimal binding conditions are adjusted and the genomic DNA is bound to the magnetic particles.

Removing residual contaminants

Contaminants are efficiently removed using **Wash Buffer I** and **Wash Buffer II**, while the DNA remains bound to the magnetic beads.

Elution

The DNA is finally eluted in **Elution Buffer**. The eluted DNA is ready-to-use in different subsequent downstream applications like:

- PCR^{*)}, RAPD, AFLP analysis
- microsatellite analysis
- genotyping
- enzymatic restriction digestion

Yield and quality of genomic DNA

The amount of purified DNA derived by the **InviMag® Plant DNA Mini Kit /KF96 w/o plastic** procedure from plant materials depends on the sample source, transport conditions, storage and sample age.

The overall yield and quality of the isolated genomic DNA is suitable for any detection system.

^{*)} The PCR process is covered by US Patents 4,683,195, and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG.

Important points before starting a protocol

Immediately upon receipt of the product, inspect the product and its components as well as the package for any apparent damages, correct quantities and quality. If there are any unconformities you have to notify Invitek Molecular in writing with immediate effect upon inspection thereof. If buffer bottles are damaged, contact the Invitek Molecular Technical Services or your local distributor. In case of liquid spillage, refer to "Safety Information" (see page 6). Do not use damaged kit components, since their use may lead to poor kit performance.

- Change pipet tips between liquid transfers. To avoid cross-contamination, we recommend the use of aerosol-barrier pipet tips.
- When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles.
- Discard contaminated gloves immediately.
- Do not combine components of different kits, unless the lot numbers are identical.
- Avoid microbial contamination of the kit reagents.
- This kit should only be used by trained personnel.

Preparing reagents and buffers

Before starting a run, bring all reagents to room temperature. Gently mix and redissolve any precipitates by warming up to 30°C. Swirl gently to avoid foaming.

Lysis Buffer P and **Elution Buffer** are ready-to-use.

Add the required amount of ddH₂O to the reaction tube containing the **Proteinase K**. Vortex for 5 s.

5x 96 DNA-extractions

Add 84 ml 99.7% Isopropanol to the **Binding Buffer A**. Mix by intensive shaking by inverting for 1 min. Shortly before use mix by inverting several times.

Resuspend each tube **Proteinase K** in 1.1 ml RNase free water, mix thoroughly until completely

Add 80 ml of 96-100% ethanol to the bottle **Wash Buffer I**

Add 140ml of 96-100% ethanol to the bottle **Wash Buffer II**

Mix thoroughly and always keep the bottle firmly closed

Reagents and equipment to be supplied by user

- Measuring cylinder (250 ml)
- Pipette and pipette tips
- Disposable gloves
- ddH₂O
- Vortexer
- 96-100% ethanol
- Isopropanol

*The **InviMag® Plant DNA Mini Kit/ KF96 w/o plastic** is validated with 2-Propanol; Rotipuran >99.7%, p.a., ACS, ISO (Order no. 6752) from **Carl Roth**.


* **Possible suppliers for Isopropanol:**

Carl Roth
2-Propanol
Rotipuran >99.7%, p.a., ACS, ISO
Order no. 6752

Applichem
2-Propanol für die Molekularbiologie
Order Nr. A3928

Sigma
2-Propanol
Order no. 59304-1L-F

Scheme of the InviMag® Plant DNA Mini Kit/ KF96 w/o plastic



Lysis Procedures

Protocol: Isolation of genomic DNA from up to 100 mg of plant material

Please read the instructions carefully and conduct the prepared procedure.

Homogenization of the starting material

Homogenize up to 100 mg of plant material by use of a pestle and mortar in combination with liquid nitrogen. Commercially available equipment for homogenization (bead mill, etc.) can be used too.

Transfer the homogenized starting material either into 1.5 ml reaction tubes (not provided) or use a 96 DWP (not provided) and add 400 µl of Lysis Buffer P and 20 µl of Proteinase K to each sample. Incubate the mixtures at 60°C for 30 min while continuously shaking.

During lysis, prefill all plates with the required buffers and appropriate volumes (see “Starting a run”, page 11). After lysis, set up the Prefilter plate I onto the Binding Plate (2.0 ml Deep Well Plate) and transfer the lysis mixtures into the filter plate. Incubate for 2-3 min at room temperature or until all lysates have passed the filter plate. The lysed samples have now been transferred into the Binding Plate.

Remove the filter plate and start the run on the KF instrument (see below).

Important: *The kit will also co-purify RNA beside DNA. For the elimination of RNA (if required), add 20 µl RNase A (10 mg/ml) to the Lysis Buffer P prior before starting the lysis procedure.*

Starting a Run on a KF96 / KFflex96 instrument

Attention: *Please be aware, that you have to prepare the **Binding Buffer A** – see instruction page: 9*

Note: *Before starting, the purification process with the KingFisher instrument please carefully read the manufacturer's manual! Resuspend/Vortex the magnetic particles (SNAP Solution) thoroughly before use!*

1. Switch on the KingFisher instrument
2. Prefill all required plates as described below:

Tip Plate: Place the KF96 Tip Comb for DW magnets on a Tip Plate (Use one provided Elution Plate as Tip Plate. These are identical.)

Binding Plate: 200 µl **Binding Buffer A**, 20 µl **SNAP Solution** to a 2 ml Deep Well Plate

Washing_Plate_1: Add 800 µl **Wash Buffer I** to a 2 ml Deep Well Plate

Washing_Plate_2: Add 800 µl **Wash Buffer II** to a 2 ml Deep Well Plate

Washing_Plate_3: Add 800 µl **Wash Buffer II** to a 2 ml Deep Well Plate

Elution Plate: Add 100 µl **Elution Buffer** to the 200 µl Elution Plate

3. Choose either the assay file “**InviMag_Plant_KF96**” (KF96 assay file) or “**InviMag_Plant_KFflex96**” (KFflex96 assay file) on the display of the corresponding KingFisher instrument and press the “START” button.
4. Insert the prefilled plates onto the right positions of the KingFisher instrument by following the specifications shown on the instrument display and confirm every loading step with the “START” button. When all prefilled plates are loaded press the “START” button to initialize the assay file. The assay file will start at the Binding Step. From this point, the instrument will continue with the purification process without any further user interaction.

The following extraction steps run automatically on the KingFisher™ System!

1. Binding of the DNA

Automatically sample mixing for 5 min. SNAPs separation. Transfer of the SNAPs to Washing Plate 1.

2. First Washing

Automatically sample mixing for 1.5 min. SNAPs separation. Transfer of the SNAPs to Washing Plate 2.

3. Second Washing

Automatically sample mixing for 1 min. SNAPs separation. Transfer of the SNAPs to Washing Plate 3.

4. Third Washing

Automatically sample mixing for 1 min. SNAPs separation.

5. Drying

Drying of SNAPs outside Washing Plate 3 for 5 minutes. Transfer of the SNAPs to the Elution Plate.

6. Elution of the DNA

Incubation of SNAPs for 10 minutes at 65°C while continuously mixing. SNAPs separation. Removal of SNAPs into Washing Plate 3 (disposal).

Important Notes:

- 1. After finishing the extraction protocol, the Elution Plate contains the extracted DNA. Store the DNA at adequate conditions. For long-term storage, we recommend to store the DNA at -20°C.*
- 2. If the extracted DNA contains carryover from magnetic particles, transfer the DNA into a new 1.5 ml reaction tube and centrifuge at maximum speed for 1 minute. Transfer the clear supernatant (containing the DNA) into a new tube.*

The eluted DNA is ready-to-use in different downstream applications. The eluted DNA can be stored for several weeks at 4-8°C or stored at -20°C for long-term storage.

For self-programming of the KF96 / KFflex96 instrument

Reagent info

Tip Plate		KingFisher 96 KF plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
-	-	-	-	
Binding Plate		Microtiter DW 96 plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Crushed sample material in Lysis Buffer P	420	-	Sample	
Binding Buffer A	200	-	Reagent	
SNAP Solution	20	-	Reagent	
Washing Plate 1		Microtiter DW 96 plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Wash Buffer I	800	-	Reagent	
Washing Plate 2		Microtiter DW 96 plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Wash Buffer II	800	-	Reagent	
Washing Plate 3		Microtiter DW 96 plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Wash Buffer II	800	-	Reagent	
Elution Plate		KingFisher 96 KF plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Elution Buffer	100	-	Reagent	

Dispensed reagents

The protocol does not contain dispensed reagents

Steps data

	Tip 1	96 DW tip comb	
	Pick-Up	Tip Plate	
	Binding Step	Binding Plate	
	Beginning of step	Precollect	No
		Release time, speed	00:00:15, Fast
	Mixing / heating:	Mixing time, speed	00:05:00, Medium
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	10
	Washing Step 1	Washing Plate 1	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Fast
	Mixing / heating:	Mixing time, speed	00:01:30, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	5
	Washing Step 2	Washing Plate 2	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Fast
	Mixing / heating:	Mixing time, speed	00:01:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	5
	Washing Step 3	Washing Plate 3	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Fast
	Mixing / heating:	Mixing time, speed	00:01:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	5
	Drying	Washing Plate 3	
		Dry time	00:05:00
		Tip position	Outside well / tube
	Elution	Elution Plate	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Medium
	Mixing / heating:	Mixing time, speed	00:10:00, Slow
		Heating temperature [°C]	65
		Preheat	Yes
	End of step	Postmix	No
		Collect count	4
		Collect time [s]	10
	Bead Removal	Washing Plate 3	
		Release time, speed	00:00:30, Fast
	Leave	Tip Plate	

Troubleshooting

Problem	Probable cause	Comments and suggestions
low amount of extracted DNA	insufficient lysis	increase lyses time, but prevent too long lyses time because this also decreases the yield Reduce amount of starting material
	incomplete elution	increase volume of Elution Buffer / increase time of elution step
	low amount of SNAP Solution	mix SNAP Solution thoroughly before addition
low concentration of extracted DNA	too much Elution Buffer	elute the DNA with a lower volume of Elution Buffer (do nt use less than 100 µl)
	incorrect storage of starting material	ensure that storage of starting material is correct avoid repeated freezing and thawing cycles of the material
degraded / sheared DNA	incorrect storage of starting material	ensure that the storage of starting material was correct avoid multiple freezing and thawing cycles of the sample material
	old material	ensure that the starting material is fresh or stored at appropriate conditions (long-term storage at -20°C) old material often contains degraded DNA
DNA does not perform well in downstream-applications (e.g. real-time PCR or PCR)	ethanol carryover during elution	increase drying time for evaporation of ethanol
	salt carryover during elution	check the Wash Buffers for salt precipitates. If any precipitates are visible, solve them by carefully warming up to 30°C ensure that the Wash Buffers are equilibrated at room temperature
low A₂₆₀:A₂₈₀ ratio from UV measurement, eluted DNA is brown colored	small part of the magnetic particles are left in the elution	centrifuge at full speed for 1 min and transfer supernatant to a new tube

Appendix

KingFisher™ BindIt Software 3.2 or higher versions

BindIt software 3.2 or higher versions were and may be used to create assay files for the KFmL, KF96/KFflex96 or KF-Duo instruments. The provided assay file(s) can either be transferred onto the corresponding workstation(s) or be started directly from within the BindIt software after assay import. Please keep in mind, that assay(s) run from within the BindIt software are not stored in the workstation memory.

Important: *Be advised that BindIt SW 3.2 or higher versions use a new unique file extension. Therefore, it is not possible to import assay files created with BindIt 3.2 or higher versions into older BindIt software versions! Please ask your local Thermo Scientific distributor for a software update.*

Note: *When creating assay files for usage with KingFisher™ instruments in combination with Microtiter Deep Well plates (e.g. Thermo Electron), it is essential to use the KingFisher™ software 3.2 or higher versions for assay development because this software version includes the correct adjustments for the microtiter plate. It is highly recommended to use Thermo Microtiter Deep Well plates with KF96 / KFflex96 / KF-Duo workstations to ensure the best purification result.*

Minimum system requirements for BindIt Software 3.2 or higher versions

PC requirements	
Supported operating systems	MS Windows XP Pro with SP3, Windows Vista SP2, Windows 7
Disk space	500 MB free disk space
Processor	Intel Pentium ≥ 1 GHz
Memory	1 GB RAM
Serial ports available	1 (for KFmL connection)
USB ports available	1 (for KF96 / KFflex96 / KFDuo connection)
Pointing device	Mouse or equivalent is required
CD-ROM drive	1
Monitor / color settings	XVGA monitor with at least 1024x768 resolution and at least a 16-bit color environment

If the actual Windows Service Packs are not installed on the corresponding lab computer, they can be downloaded from the Microsoft web pages: <http://www.microsoft.com/>

General notes on handling DNA

Nature of DNA

The length and delicate physical nature of DNA requires careful handling to avoid damage due to shearing and enzymatic degradation. Other conditions that affect the integrity and stability of DNA include acidic and alkaline environments, high temperature, and UV irradiation. Careful isolation and handling of high molecular weight DNA is therefore required to ensure compatibility with various downstream applications. Damaged DNA may perform poorly in applications such as genomic Southern Blotting or long-template PCR.

Storage of DNA

A working stock of DNA can be stored at 2-8°C for several weeks. For long-term storage, DNA should be stored at -20°C, but storing at -20°C may cause shearing, particularly if the DNA is exposed to repeated freezing and thawing cycles.

Note that the solution, in which the nucleic acid is eluted in, will affect its stability during storage. Pure water lacks buffering capacity and an acidic pH may lead to acid hydrolysis. Tris or Tris-EDTA buffer contains sufficient buffering capacity to prevent acid hydrolysis.

Drying, dissolving and pipetting DNA

Avoid overdrying of genomic DNA after ethanol precipitation. We highly recommend to air dry than to use a vacuum, although vacuum drying can be used with caution.

Avoid vigorous pipetting. Pipetting of genomic DNA through small tip openings can cause shearing or nicking. One way to decrease shearing of genomic DNA is to use special tips that have wide openings designed for pipetting genomic DNA.

DNA yield

The amount of purified DNA from the plant material depends on sample source, transport conditions, storage and age of the sample.

Ordering information

InviMag® Plant DNA Mini Kit/ KF96 w/o plastic	7437300250	5 x 96 preps
Invisorb® DNA Plant HTS 96 Kit/ C	7037300300	4 x 96 preps
Invisorb® DNA Plant HTS 96 Kit/ C	7037300400	24 x 96 preps
Invisorb® Spin Plant Mini Kit	1037100300	250 purifications

Possible suppliers for Isopropanol:

Carl Roth
2-Propanol
Rotipuran >99.7%, p.a., ACS, ISO
Order no. 6752

Applichem
2-Propanol für die Molekularbiologie
Order no. A3928

Sigma
2-Propanol
Order no. 59304-1L-F

Ordering information (KingFisher™ 96 and consumables)

Cat.no	Description
5400500	KingFisher 96, magnetic particle processor, 100-240V, 50/60Hz (including one magnetic head)
24073430	KingFisher 96 head for Deep Well Plates
97002514	KingFisher 96 tip comb for a PCR magnet head / plates, 8 x 10 pcs / box
97002524	KingFisher 96 tip comb for KF magnets / plates, 10 x 10 pcs / box
97002534	KingFisher 96 tip comb for DWP magnets / DWP 10 x 10 pcs / box
97002540	KingFisher 96 KF plate (200 µl) 48 plates / box
95040450	Microtiter deep well 96 plate (2 ml), 50 plates / box

INVITEK
Molecular

Invitek Molecular GmbH
Röbert-Rössle-Str. 10
13125 Berlin

Phone: +49 30 94 89 29 01
Fax: +49 30 94 89 29 09
info@invitek-molecular.com

www.invitek-molecular.com

0067437300 V-01-2020