SalivaGene Collection Module II

DEUTSCH

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SalivaGene Collection Module II - Instructions

For In Vitro Diagnostic Use*

*) Compliance with EU Directive 98/79/EC on in vitro medical devices (Not for in vitro diagnostic use in countries where the EU Directive 98/79/EC on in vitro medical devices is not recognized.)

present and easy to do. Continue rinsing.

LABEL LEGEND:

Consult operating instructions

Store at room temperature 15 - 30°C.

Do not reuse

Expiry date

1035212200

CE Marking

In vitro diagnostic medical device

Invitek Molecular GmbH Robert-Roesse-Str. 10 13125 Berlin, Germany

+49 30 9489 2908 info@invitek-molecular.com

www.invitek-molecular.com

INTENDED USE: Collection of human saliva samples for DNA stabilization and extraction (using the SalivaGene DNA Extraction Kits). The DNA in SalivaGene stabilized samples is stable for 1 year at room temperature within the expiry date. The product is intended for use under supervision of professionals only, such as technicians, physicians and biologist trained in molecular biological techniques.

CONTENTS: Collection tube, Stabilizer tube. The stabilizer tube contains approx. 2 ml of liquid Saliva DNA Stabilizer reagent.

LIMITATIONS OF THE PROCEDURE: Do not eat, drink, smoke or chew gum at least 30 min prior to saliva collection. Collecting 2 ml of saliva may require several minutes. Putting some grains of sugar on the tongue helps to produce saliva and does not interfere with the stabilization.

1. Unscrew lid from Collection tube and put it aside for further use.

2. Rub cheeks against teeth intensely for 30 sec.

3. Collect saliva to indicated fill level, avoid making and measuring air bubbles.

4. Unscrew and discard lid from Stabilizer tube. Pour Saliva DNA Stabilizer reagent into Collection tube.

5. Screw lid tightly onto the Collection tube again.

6. Shake tube for 15 sec to mix saliva and Saliva DNA Stabilizer reagent.

STORAGE: Store at room temperature 15 - 30°C (59 - 86°F).

SHIPMENT: Ship collected saliva according to the instructions of SalivaGene Collection Module II provider. Ship according to locally applicable guidelines and regulations such as UN 3373 (Biological Substance, Category B).

SAFETY PRECAUTIONS: Keep out of the reach of children. Small cap: choking hazard. Do not use the SalivaGene Collection Module II if content of the tube is spilled. Do not drink the Saliva DNA Stabilizer reagent.

Wash with water if Saliva DNA Stabilizer reagent or Stabilizer-Saliva reagent mix comes in contact with skin or eyes.

WARNING: H319: Causes serious eye irritation. H412: Harmful to aquatic life with long lasting effects. P280: Wear protective gloves/protective clothing/ eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if

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